



## Microsoft Update and PHP Update Monthly Notification

**Summary:** This notification contains information about Microsoft security updates released this month and approved for implementation on the Gyrus ACMI, Inc., ES-IP, ES4K, LiveStream™, nCare™, and VaultStream™ platforms. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the EasySuite™ HD, ES4K, and ES-IP iOR systems, nCare MX/RX recorders, EasyLink™ routers, LiveStream StatusBoard Engine, and VaultStream systems and its KE Modules.

Also included is information on the latest PHP updates that are approved for the v10.x and up VaultStream systems.

**Products:** ES-IP and ES4K iOR systems and nCare recorders, connected to a network and VaultStream server software, as well as the LiveStream StatusBoard Engine.

**Date:** December 12, 2023

### **Critical Security Bulletins (Obtained from Microsoft):**

**NOTE:** Gyrus ACMI, Inc. **does not approve** MS updates **KB975560, KB975562, KB2631813, or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

### **Approved Updates:**

<b>December 2023 Updates</b>		
<b>KB Number</b>	<b>Platform</b>	<b>Title</b>
<b>VaultStream v10.x Updates – Windows Server 2012</b>		
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.366.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.353.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.322.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.304.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.260.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.212.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.202.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.149.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.137.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.111.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.87.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.70.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.30.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1710.0)

KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1703.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1664.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1656.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1618.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1602.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1583.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1565.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1546.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1538.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1491.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1481.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1436.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1431.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1392.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1356.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1336.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1303.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1255.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1223.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1212.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1170.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1160.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1115.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1056.0)

KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.1011.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.963.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.929.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.908.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.883.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.856.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.841.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.798.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.754.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.696.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.661.0)
KB2310138	VaultStream 2012	Security Intelligence Update for Microsoft Security Essentials - (Version 1.403.654.0)
<b>VaultStream v10.x Updates – Windows Server 2016</b>		
KB5033373	VaultStream 2016	2023-12 Cumulative Update for Windows Server 2016 for x64-based Systems
<b>VaultStream v10.4 Updates – Windows Server 2019</b>		
KB5033371	VaultStream 2019	2023-12 Cumulative Update for Windows Server 2019 (1809) for x64-based Systems
<b>nCare MX v10.2.3</b>		
KB5033373	nCare MX	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB2267602	nCare MX	Security Intelligence Update for Microsoft Defender Antivirus
KB2267602	nCare MX	Security Intelligence Update for Microsoft Defender Antivirus
<b>nCare MX v10.4</b>		
KB5033373	nCare MX	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB2267602	nCare MX	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.3887.0)
KB2267602	nCare MX	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.1963.0)
<b>ES4K v10.4 ES4K Controller (a.k.a. BlizCon)</b>		
KB5033373	ES4K Controller	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB2267602	ES4K Controller	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.3887.0)
KB2267602	ES4K Controller	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.1963.0)

<b>ES4K v10.4   ES4K Video Manager (a.k.a. nCare RX2)</b>		
KB2267602	ES4K Video Manager	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.3887.0)
KB2267602	ES4K Video Manager	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.1963.0)
<b>ES4K v10.4 LiveStream StatusBoard 2</b>		
KB5033373	LiveStream StatusBoard 2	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB2267602	LiveStream StatusBoard 2	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.3887.0)
KB2267602	LiveStream StatusBoard 2	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.1963.0)
<b>MedPresence Tower v10.4</b>		
KB5033373	MedPresence Tower	2023-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB2267602	MedPresence Tower	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.403.404.0)
KB2267602	MedPresence Tower	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.391.1963.0)
<b>ES-IP iOR system</b>		
KB5033372	ES-IP	2023-12 Cumulative Update for Windows 10 Version 21H2 for x64-based Systems
KB890830	ES-IP	Windows Malicious Software Removal Tool x64 - v5.119
<b>MAJ-2533 (UCES-4)</b>		
KB5033371	MAJ-2533	2023-12 Cumulative Update for Windows 10 Version 1809 for x64-based Systems
<b>CV Connectivity with KE v3.1 (VaultStream with Interfaces License)</b>		
KB5032189	KE v3.1	2023-11 Cumulative Update for Windows 10 Version 22H2 for x64-based Systems
KB5032339	KE v3.1	2023-11 Cumulative Update for .NET Framework 3.5, 4.8, and 4.8.1 for Windows 10 Version 22H2 for x64
KB890830	KE v3.1	Windows Malicious Software Removal Tool x64 - v5.119
KB5032197	KE v3.1	2023-11 Cumulative Update for Windows Server 2016 for x64-based Systems
KB5032391	KE v3.1	2023-11 Servicing Stack Update for Windows Server 2016 for x64-based Systems
KB5032196	KE v3.1	2023-11 Cumulative Update for Windows Server 2019 (1809) for x64-based Systems
KB5032337	KE v3.1	2023-11 Cumulative Update for .NET Framework 3.5, 4.7.2 and 4.8 for Windows Server 2019 for x64
<b>PHP Updates</b>		
V10.2.2 was Tested with PHP version 7.3.33 V10.4 was Tested with PHP version 7.4.30 7.3.33 is available upon request for v10.2.2 VaultStreams only. 7.4.33 is available upon request for v10.4 VaultStreams only.		

For further questions concerning this notice, please call Technical Support at 1-800-848-9024.

Microsoft and Windows are trademarks of the Microsoft group of companies.



**Gyrus ACMI, Inc.**  
 One 9600 Louisiana Ave. North, Brooklyn Park, MN 55445 U.S.A.  
[www.olympusamerica.com/olympus-systems-integration](http://www.olympusamerica.com/olympus-systems-integration)