



December 9, 2019

Re: Compatibility of Acecide-C with Olympus ultrasonic probes

Dear Health Care Professional,

This letter is in regards to the compatibility of Acecide[®]-C high level disinfectant with Olympus ultrasonic probes.

The following Olympus ultrasonic probes are compatible with Acecide[®]-C and can be reprocessed in the Olympus OER-Pro automated endoscope reprocessor:

**UM-2R
UM-3R
UM-S20-17S
UM-S20-20R
UM-S30-20R**

**UM-S30-25R
UM-BS20-26R
UM-G20-29R
UM-DG20-31R
UM-DP12-25R**

**UM-DP20-25R
RU-12M-R1
RU-75M-R1**

Please note that UM models require the MAJ-1764 connecting tube and RU models require the MAJ-1515 connecting tube when reprocessing within the OER-Pro.

For a list of all devices compatible with Acecide[®]-C within the OER-Pro, please refer to the Olympus document titled *List of Compatible Endoscopes/Connecting Tubes <OER-Pro>*. Inclusion of a device within this document infers compatibility with Acecide[®]-C high level disinfectant.

If you have any additional questions, please contact your local Olympus representative or the Olympus Technical Assistance Center at 1-800-848-9024 (United States) or 1-800-387-0437 (Canada).

Sincerely,

Olympus